



Vancouver Diagnostic Services Immunochemistry Referral Testing Services

Direct Antiglobulin Test (DAT)

BC_REF-03

TEST DESCRIPTION

This test is performed to determine whether a patient's red blood cells have been sensitized in vivo with immunoglobulin, complement, or both. The DAT is used most commonly to investigate possible hemolytic transfusion reactions, Hemolytic Disease of the Fetus and Newborn (HDFN), autoimmune hemolytic anemia, and drug-induced immune hemolysis. The clinical significance of a DAT result should take into consideration the patient's clinical history, diagnoses, and other laboratory test results.

SPECIMEN AND REQUISITION REQUIREMENTS

Specimen(s)

- Minimum of three (3) 6-7ml (lavender) tubes, unless pediatric patient, mixed thoroughly by gentle agitation.
- Label specimen with the required minimum information: patient's last name, first name, PHN or unique identifier and date of collection.

Complete Requisition (must include)

- Patient's last name, first name, date of birth and PHN or unique identifier
- Facility name, complete address, phone and fax number
- Phlebotomist ID information
- Date of collection

Requisition(s)

- *Antibody Investigation Requisition_BC*

PRE-SHIPING STORAGE

Recommended Refrigeration 1-10°C.

SHIPPING INSTRUCTIONS

Submit samples as soon as possible after collection.

Shipping

- Ship in a container that will maintain temperature at $\geq 1^{\circ}\text{C}$.
- Select shipping method for container to arrive at testing site within 48 hours.

Note: *Protect from freezing.*

SEND TO

Canadian Blood Services
BC & Yukon Centre
Diagnostic Services Laboratory
4750 Oak Street
Vancouver, BC V6H 2N9
Tel: 604-707-3434
Fax: 604-874-6582